



VetX DR3326M

Full HD XXL Imaging

The mobile X-Ray DR-System is fitted into a sturdy case with ergonomically integrated electronics. The DR-plate is protected separately for maximum mobility and optimal use of space. The mobile X-Ray DR-System is an extremely light-weight yet full-scale DR-System.

It consists of a virtually indestructible case, the digital imaging plate, a powerful business-class notebook and the entire electronic-set incl. cables. Consequently this system is easily carried or comfortably pulled along on its wheels.

The system acquires your x-ray images in a cycle of approx. 4 seconds and transfers them to your display in less than approx. 2 seconds. The silicium plate-technology is less susceptible to temperature, humidity and shock, making it the perfect DR-system for equine application.



Winden 26, 85304 Ilmmünster
Germany

eMail info@xray-is.de fon +49 8441 / 479 229-3 (Sales)
support@xray-is.de +49 8441 / 479 229-2 (Support)

internet www.xray-is.de fax +49 8441 / 479 229-9

SYSTEM BENEFITS

- Panel 11 x 13' (26 x 33 cm)
- Sturdy and designed for veterinarians
- Highest dynamic (16 bit raw data)
- Low susceptibility to faults
- Small weight, panel ca. 3.4 kg only
- High scan- and processing-speed

CONVINCING TECHNOLOGY

- Sand-, water- and air-tight transport-case
- Made from HPX® - Virtually indestructible
- Soft-Grip handle with press-pull-latches
- Sturdy in-case wheels, telescopic handle

Diagnostic Software

The ease of use and, above all, the perfect adjustment to daily routines in a veterinary practice of our software is well-known.

The software has been designed for the specific needs of veterinary requirements, and supports your every day image processing and findings perfectly. Our image-filters have been fine-tuned and optimized for the veterinary market; moreover, on installation, they are adjusted to your specific wishes and your x-ray generator in particular, of course. There is a difference to others - flexibility helps us, and You, to keep the edge!

The Diagnostic Workstation has been sold 9000 times in ca. 67 countries for some reason.



Your dealer:

SOFTWARE HIGHLIGHTS

Input screen specifically designed for veterinary application - GDT-interface to PMS-systems (e.g. VetInf) - direct image acquisition from various modalities such as ultrasound, endoscopy, CT, etc - image processing filters specifically designed for veterinary application - special gray-scale look-up-tables (LUTs) for equine x-ray imaging - measurement-tools e.g. hip dysplasia (dogs) - CT-Series 3D-Imaging - PACS-Interface and DICOM-Print, DICOM-Email, DICOM-Worklist, etc - additional network-workstations (viewer) - individual filter adjustments (via sliders) Currently the VetX Diagnostic Workstation is available in the following languages: German, English, French, Italian, Polish, Russian, Spanish and Chinese

SYSTEM FEATURES

- DR-Panel 11 x 13' (26 x 33 cm)
- 17.3" Business-Notebook 500GB/4GB, USB2.0, Win 7/32Bit (features always according to current standard)
- VetX Diagnostic Workstation Standard, Patient-CD Interface (connect) to DR System incl. special Vet-Filter
- Installation & configuration, Hardware-briefing
- Training on VetX Diagnostic Workstation (ca. 2 hrs.)
- System in HPX-case, plate in separate padded case
- incl. trigger-box for sync. of x-ray generator and DR-panel
- 1 year warranty on panel or according to manufacturer's warranty!
- Dimensions (W/D/H)/Weight:
System package 56x35x23 cm/ca. 18kg

PANEL TECHNICAL DATA

Pixel size	127µm
Active pixel area	2560 x 2080mm
AD converter	14 bit , 16384 gray scales
Resolution	3.9 lp/mm
HR-image display	ca. 2 seconds

OPTIONAL

- Wireless trigger-box for sync. of x-ray generator and DR-Plate
- Wireless data transfer with beltbag for accu