

X-Ray Generators

X-Ray Generators for Out-Patient Application!

- Portable and powerful
- Various models available
- Dual laser-pointer for easy distance-measuring
- All 'Orange' generators equipped with sturdy transport-case

ORANGE

Orange 8016HF/9020HF Portable



Orange 1025HF/1040HF/1060HF Portable



ORANGE PORTABLE

High Frequency Generator
Fixed Anode
SID Scale length: Build-in with DUAL Laser beam

Orange 801 6HF Portable

80kV / 16mA, (Superior tube), Focal Spot: 1.0 mm
kV Range: 40 - 80kV, mAs Range: 0.32 - 30 mAs
Weight: 6.2 kg, Max. Output: 1.2Kw @ 80 kV

Orange 9020HF Portable

90kV / 16mA, Focal Spot: 0,8 mm
kV Range: 40 - 90kV (1kV step), mAs Range: 0.32 - 30 mAs
Weight: 6.9 kg, Max. Output: 1.35 Kw @ 80 kV

Orange 1025HF Portable

100kV / 25mA, Focal Spot: 1,2 mm
kV Range: 40 - 100kV (1kV step), mAs Range: 0.32 - 50 mAs
Weight: 9.9 kg, Max. Output: 1.6 Kw @ 80 kV

Orange 1040HF Portable

100kV / 40mA, Focal Spot: 1,2 mm
kV Range: 40 - 100kV (1kV step), mAs Range: 0.32 - 100 mAs
Weight: 11,5 kg, Max. Output: 2,4 Kw @ 80 kV

Orange 1060HF Portable

100kV / 60mA, Focal Spot: 2,0 mm
kV Range: 40 - 100kV (1kV step), mAs Range: 0.32 - 100 mAs
Weight: 12,9 kg, Max. Output: 3,2 Kw @ 80 kV